


Clinic | Care | Homecare | Assist







# Design details



Reverse-Trendelenburg position



Adjustment range from 30 - 80 cm



Electronic switch mode power supply unit integrated immediately after the mains plug



Handset with selective locking device

# Allura II -Innovative care beds, including for heavier occupants

## Advantages

## **Quality and durability**

The "Allura II" care bed is extremely sturdily constructed to stand up to the rigours of daily care. High-quality components make this a particularly durable, sturdy and low maintenance care bed that significantly reduces follow-on costs.

## Low entry height and safe working load of 250 ka

The Allura II care bed series is available with mattress base dimensions of 90 x 200, 100 x 200 or 120 x 200 cm. With a safe working load of 250 kg, the care bed is also ideally suited to heavier persons. The mattress base height can be adjusted from 30 - 80 cm. The exceptionally large adjustment range of 50 cm ensures a high degree of nursing comfort.

## 24-volt drive system as standard

By using the 24-volt drive system, the prescribed intervals for leakage current testing are extended to such a degree that leakage current testing is no longer required during or beyond the expected service life of the bed. Please refer to the instruction manual for the length of this period. The first leakage current test has already been conducted by the manufacturer. The 24-volt extra-low voltage saves electricity and enhances the safety of occupants.

## Safety side protective height of 38 cm

Mattresses and anti-decubitus systems up to a maximum thickness of 16 cm can be used in the care bed without the need for adaptable safety side height extensions. Expensive height extensions are unnecessary in most cases.



Compact delivery and storage on storage aid/ transportation aid (similar to illustration)

# Allura II care bed

# Features

- Mattress base: 90 x 200 cm (Allura II), 100 x 200 cm (Allura II 100) or 120 x 200 cm (Allura II 120), 4-section, mattress base height adjustable from 30 - 80 cm, backrest, thigh rest and mattress base height are adjustable with electric motors
- Reverse-Trendelenburg position can be set by tilting the mattress base
- Metal parts finished with polyester powder coating in the colour topaz
- Headboard and footboard panels in beech decor
- Handset with selective locking function
- Safe working load: 250 kg
- 24-volt drive system

Manufactured to: standard DIN EN 60601-2-52: 2010

# Accessories

(available at extra cost)



Bed extension (also for retrofitting)



Reading lamps



Foam leather cover



Tray for Allura II and Allura II 100



Wooden panels for the chassis (also for retrofitting)

# Order numbers and external dimensions

51.5024.33	Allura II care bed (wooden safety sides), length 224 cm, width 101 cm

- 51.5524.33 Allura II 100 care bed (wooden safety sides), length 224 cm, width 111 cm
- 51.5724.33 Allura II 120 care bed (wooden safety sides), length 224 cm, width 131 cm
- 19.1400.1 Lifting pole (included in supply package)

 Burmeier GmbH & Co. KG

 Pivitsheider Straße 270, 32791 Lage, Germany

 Phone:
 +49 (0) 5232 9841 - 0

 Fax:
 +49 (0) 5232 9841 - 41

 Email:
 info@burmeier.com

 www.burmeier.com



# CE