

Inovia II

Clinic | Care | **Homecare** | Assist




The care bed with
innovative ideas



Practical solid wood grip rail strikes a home-like note

Protective height of safety sides of 41 cm for mattresses up to 19 cm

Sliding backrest



Low wooden footboard is an attractive feature

Castors can be locked in pairs in every mattress base position

Easy-clean mattress base slats

Inovia II - The care bed with innovative ideas

Advantages

Mattress base height adjustment from 22 - 77 cm

The innovative 55 cm adjustment range makes this care bed ideal as a low-entry bed for fall protection.

24-volt drive system as standard

By using the 24-volt drive system, the prescribed intervals for leakage current testing are extended to such a degree that leakage current testing is no longer required during or beyond the expected service life of the bed. Please refer to the instruction manual for the length of this period. The first leakage current test has already been conducted by the manufacturer. The 24-volt extra-low voltage saves electricity and enhances the safety of occupants.

Safe working load of 225 kg

Also suitable for heavier persons.

Safety side protective height 41 cm

Mattresses and anti-decubitus systems up to a maximum thickness of 19 cm can be used in the care bed without the need for adjustable safety side height extensions. Expensive height extensions are unnecessary in most cases.

Mattress base with sliding backrest in line with German Nursing Association recommendations

The retracting backrest protects the resident against uncomfortable shearing forces while the mattress base is adjusted and thus contributes to the prevention of pressure ulcers. The retraction of the backrest enlarges the seating area for the resident and increases his comfort significantly.

Design details



Solid wood grip rail



Castors can be locked in pairs in every mattress base position



Handset with selective locking



Sliding backrest



Mattress base adjustable from 22 - 77 cm, adjustment range = 55 cm

Inovia II care bed

Features

- Mattress base 90 x 200 cm or 100 x 200 cm, 4-section, with sturdy, easy-clean steel slats, mattress base height adjustable from 22 - 77 cm, backrest, thigh rest and mattress base height are adjustable with electric motors. Reverse-Trendelenburg position adjustable by tilting the mattress base
- Mattress base with sliding backrest in line with German Nursing Association recommendations
- to the prevention of pressure ulcers. The receding backrest enlarges the seating area for the resident and increases his comfort significantly.
- Castors can be locked in pairs in every mattress base position
- Safe working load: 225 kg
- Wooden surround in beech decor/solid beech
- Footboard with solid wood grip rail (also headboard as an option)
- Integrated safety sides, protective height 41 cm for mattress thicknesses of up to 19 cm
- Innovative, home-like appearance, with low wooden footboard despite 55 cm height adjustment range
- 24-volt drive system



Manufactured to: standard IEC 60601-2-52:2009

Subject to technical changes without notice.
Last updated: 03/2017, Version 5, Item no. 600142.

Accessories

(available at extra cost)



Tray



Foam leather cover
for safety sides



Reading lamps

- Bed extension
- Wall distance holders
- Wooden surround in Real Cherry or Royal Maple

Order numbers and external dimensions

- 50.5522.24 **Inovia II** care bed, length 220 cm, width 105 cm
- 55.5522.34 **Inovia II 100** care bed, length 220 cm, width 115 cm
- 19.1400.1 Lifting pole

Burmeier GmbH & Co. KG

Pivitsheider Straße 270, 32791 Lage, Germany

Phone: +49 (0) 5232 9841 - 0

Fax: +49 (0) 5232 9841 - 41

Email: info@burmeier.com

www.burmeier.com

