

Westfalia-Care II

Clinic | Care | **Homecare** | Assist



The versatile
care bed

Westfalia-Care II - The versatile care bed

Manually adjustable
lower leg rest

Integrated bed extension
and lower leg extension

Durable, low maintenance
LINAK drive system
(OpenBus system)

225 kg safe
working load

Central
castor locking

Integrated
linen holder

Lowest mattress base position
of only about 26 cm

Sliding
backrest

Advantages

Mattress base height adjustment from 26 - 89 cm

The innovative 63 cm adjustment range makes this care bed ideal as both a low-entry bed and a standard care bed. The highest mattress base position of 89 cm allows care staff to work at a comfortable, back-friendly height. The need to define a single use is now a thing of the past – this bed is a real multi-talent.

Safe working load of 225 kg

Also suitable for heavier people.

Design

The innovative and timeless design emphasises the Burmeier philosophy, with its pleasing appearance, uncompromising high-quality and durability.

24-volt drive system as standard

By using the 24-volt drive system, the prescribed intervals for leakage current testing are extended to such a degree that leakage current testing is no longer required during or beyond the expected service life of the bed. Please refer to the instruction manual for more information about this. The first leakage current test has already been conducted by the manufacturer. The 24-volt extra-low voltage saves electricity and enhances the safety of occupants.

Mattress base with sliding backrest in line with German Nursing Association recommendations

The retracting backrest protects the resident against uncomfortable shearing forces while the mattress base is adjusted and thus contributes to the prevention of pressure ulcers. The retraction of the backrest enlarges the seating area for the resident and increases his comfort significantly.

Design details



Central castor
locking possible in
any mattress base
position, castors
concealed
with cover



Integrated linen
holder (only
available with
the Primus
footboard shown)



Integrated
bed extension
and lower leg
extension



Electronic switch
mode power supply
unit integrated
straight after the
mains plug

Westfalia-Care II care bed

Features

- Mattress base 90 x 200 cm, one-piece (not dismantlable), with sturdy, easy-clean steel slats, in the colour titanium metallic
- Mattress base height adjustment range 26 - 89 cm, backrest, thigh rest and bed height adjustable with electric motors
- Concealed castors, can be centrally locked in every mattress base position
- Mattress base with sliding backrest in line with German Nursing Association recommendations
- Safe working load: 225 kg
- Wooden surround in beech decor/solid beech
- Footboard with integrated bed extension (20 cm), wooden surround with integrated safety sides
- 24-volt drive system

Manufactured to: standard IEC 60601-2-52:2009

Subject to technical changes without notice.
Last updated: 03/2017, Version 4, Item no. 600425.



Accessories

(available at extra cost)



Tray



Wall deflection roller with holder



Reading lamps




Wooden head and footboards:
Standard (as standard), Comodo
with round posts (optional),
Prinzino with handle bar
(optional)



Wooden head and footboards:
Primus headboard closed/
footboard with handle bar,
with or without integrated
linen holder



Foam leather cover



Longer safety sides for the
integrated bed extension

Order numbers and external dimensions

50.5432.12 **Westfalia-Care II** care bed, length 212 cm, width 102 cm

19.1400.1 Lifting pole

Burmeier GmbH & Co. KG

Pivitsheider Straße 270, 32791 Lage, Germany

Phone: +49 (0) 5232 9841 - 0

Fax: +49 (0) 5232 9841 - 41

Email: info@burmeier.com

www.burmeier.com

